

CANION USER GUIDE

(V3.3, Updated February 2015)

CANION is program for analyzing the distribution of ions, solute atoms, or water molecules surrounding (or belonging to) a nucleic acid. The data for CANION can come in three formats: (1) an Amber format molecular dynamics trajectory that has already been analyzed with Curves+ (*file.cdi*); (2) a density matrix in Gaussian cube format (*file.cub*); (3) a simplified snapshot format with a specified number of ion/atom Cartesian coordinates per snapshot. (*file.pts*) We will start by describing how to generate data for Canion using Curves+ and the options for analyzing this data in Canion. We will then describe the differences when using the other input options.

In the following text we will refer only to analyzing "ions", but the same analysis can equally be applied to atoms other than ions, including the oxygen atoms of solvent water molecules, or atoms belonging to the nucleic acid solute molecule.

First step: CURVES+ analysis options

The ions to be analyzed in CURVES+ are specified in the file `standard_i.lib`. Each line contains a name (in single quotes) and a formal charge. Lines starting with # are treated as comments. Here is an example:

```
# Ion library
# Each line gives the ion/atom name (in single quotes) and its
# formal charge
# Up to 40 ion/atom types are allowed in this library file
'K+' 1
'Cl-' -1
'P' -1
'C1*' 0
'O' 0
```

Note sugar ring atoms names in Amber input containing quotes are changed to stars by Curves+ (e.g. `C1'` → `C1*`). Therefore, only starred sugar atom names are used in `standard_i.lib`.

To analyze ion positions, use the Curves+ namelist option `ions=.t.` (see Curves+ user guide). For each snapshot composing a trajectory, Curves+ calculates a helical axis and then determines the position of the ions with respect to the axis. This involves finding the point on the axis that corresponds to the shortest perpendicular distance to the ion. The ion position is then defined by the distance of this point along the axis in units of base pair steps (**d**), the length of the perpendicular vector (**r**) and the angle of this vector with respect to the vector of the local axis frame pointing towards the 5' → 3' strand (**a**). This choice places the center of the minor groove in B-DNA at roughly at **a** = 90° and the center of the major groove roughly at **a** = 270°.

The ion position are output by Curves+ in .cdi file containing the curvilinear helicoidal coordinates (CHC) of each ion in each snapshot.

Note that ions will only be analyzed within a radius of 30 Å from the helical axis of the nucleic acid. Beyond this distance the solute molecule has no significant influence of the bulk ion distribution. Ions lying beyond the ends of the helical axis are similarly ignored (a necessary choice given the CHC used to analyze ion distributions)

Analyzing the CHC data in Canion also requires an average helical axis. This can be obtained by generating an average nucleic acid structure for the trajectory (for example, using ptraj) and analyzing this structure in Curves+ using the option axfrm=.t. (see Curves+ user guide). This will generate an .afr file containing helical axis frames at each base level of the nucleic acid. Alternatively, Canion can generate a uniform linear helical axis internally, using helical rise and twist values set by the user (see below).

CANION analysis using data from Curves+

CANION reads the ion CHC data from the .cdi file generated by Curves+. It optionally reads an average helical axis (generated as a set of frames by Curves+ as described above).

CANION namelist variables

| Name | Type | Default | Meaning |
|--------|-----------|--------------|---|
| lis | character | | O/P .lis list file |
| dat | character | | I/P file with extension .cdi, .cub or .pts |
| axfrm | character | | I/P .afr file name |
| prop | character | | read prop .ser data for filtering |
| type | character | * (i.e. all) | ion type to be analyzed |
| seq | character | * (i.e. all) | base sequence to be analyzed |
| seqin | character | | For cube/points I/P, give 1st strand sequence |
| solute | character | | I/P .DNA pdb file for calculating solute volume |
| sris | real | 3.38 Å | standard rise if no I/P .afr file |
| stwi | real | 34.5° | standard twist if no I/P .afr file |
| grid | real | 1.0 Å | Grid spacing for 3D distribution |
| alow | real | 0.0° | Lower angle sampling limit |
| ahig | real | 360.0° | Upper angle sampling limit |
| dlow | real | 1.0 bp | Lower bp sampling limit (1 → nlev) |
| dhigh | real | 500.0 bp | Upper bp sampling limit (1 → nlev) |
| rlow | real | 0.0 Å | Lower radius sampling limit |
| rhig | real | 30.0 Å | Upper radius sampling limit |
| pmin | real | -500 | minimum filtering limit using I/P .ser data |
| pmax | real | 500 | maximum filtering limit using I/P .ser data |

| Name | Type | Default | Meaning |
|---------------|---------|---------|---|
| iprop | integer | 0 | base pair level property to be used with pmin and pmax for filtering (1 → nlev) |
| itst | integer | 0 | First snapshot to analyze |
| itnd | integer | 0 | Last snapshot to analyze |
| itdel | integer | 1 | Spacing of snapshots to analyze |
| istep | integer | 0 | If >0 O/P total number of ions every isteps to lis.stp file |
| circ | logical | .f. | Set true for minicircle analysis |
| rmsf | logical | .f. | Set true for rmsf calculation per "ion" |
| series | logical | .f. | Set true for series output |

Very important note: (1) namelist input starts with &inp and ends with &end; (2) namelist data lines must not start on column 1; (2) it is not necessary to put quotes around character variables.

Further explanations

1) Selecting the ions to analyze:

type selects the ions (or atoms) to be analyzed. It can be a name (corresponding to one of the entries in the standard_i.lib file) used during the Curves+ analysis, or a number (1, -1 or 0 zero - meaning analyze all cations, anions or neutral species), or "*" meaning analyze all available data.

2) Selecting the sequence elements to analyze:

seq selects base sequences to analyze, it contains a character string made up of A, G, C, T, R, Y, S, K, W, M or *. (Note: R=A/G, Y=C/T, S=G/C, K=G/T, W=A/T, M=A/C). The I/P sequence will be searched for any occurrences of this string. If the string has an odd number of characters the analysis is limited to the central base pair of the string (extending ½ a base pair step in each direction). If the string has an even number of characters the central base pair step is analyzed. **seq=*** will analyze the whole I/P oligomer. Note that the **seq** string is searched for in both strands. If it occurs in the second strand the 2D results are dyad inverted so they are compatible with first strand results. Which axis segments of the nucleic acid finally are analyzed are indicated in the Canion .lis file. The example below shows a fragment of the Canion output for an sample oligomer using seq=TT:

```

Preselect = C      T      T      C      T      A      T      A      A
             |-----|^^^|-----|-----|-----|-----|vvvvvv|vvv
             A      A      G      G      C      T      G
             vvv|vvvvvv|-----|-----|-----|-----|

```

- for this oligomer "TT" selects four base pair steps, one in the 1st strand ('^') and three in the 2nd strand 'v'). Each step is divided into six subdivisions.

3) Selecting the oligomer fragment to analyze:

dlow, **dhig** can also be used to select a segment of the oligomer, between base pair levels 1 and nlev (where nlev is the number of base pairs in the oligomer). With the default values (**dlow**=**dhig**=0) the entire oligomer is treated. **dlow** and **dhig** (where dlow < dhig) can be set to any real number between 1 and nlev. Note that **seq** and **dlow/dhig** can be used together, but be careful to check the Canion output file to see what has actually been selected.

4) Selecting the radial and angular ranges to analyze:

rlow and **rhig** fix the minimal and maximal distances from the helical axis taken into account in the analysis. It may be useful to note that in canonical B-DNA the phosphorus atoms are roughly 10.25 Å from the axis and rhig=10.25 is generally used to limit the analysis to the grooves of the double helix.

alow and **ahig** fix the minimal and maximal angle range to analyze. Note that angle around the helical axis are measured from 0° to 360°, where 0° corresponds to the 5'-3' strand direction in canonical B-DNA and 90° to the center of the minor groove. In canonical B-DNA, the phosphorus atoms lie at 33° and 147°. If **alow** ≤ **ahig**, angles between these values are accepted, i.e. **alow**=33° and **ahig**=147° defines the minor groove. If **alow** > **ahig**, angles outside the **ahig** to **alow** range are accepted, i.e. **alow**=147° and **ahig**=33° defines the major groove (and accepts angles 0°-33° and 147°-360°).

5) Selecting the snapshots to analyze:

itst, **itnd**, **itdel** enable the I/P data to be filtered in order to study part of the trajectory or to increase the spacing between the snapshots analyzed. **itst** sets the first snapshot analyzed, **itnd** sets the last snapshot and **itdel** sets the spacing between analyzed snapshots. Setting **itst**=n (n > 0) and leaving **itnd**=0, will analyze a snapshot number n.

6) Filtering the analysis on the basis of existing data:

The namelist options **prop**, **pmin**, **pmax** trigger reading a **prop.ser** file from a Canion analysis (which should contain 1 line of data per snapshot, with one data value per base pair level in the oligomer being analyzed). **iprop** determines the base pair level to be used for filtering. If, for a given snapshot, the corresponding data value falls outside the pmin/pmax limits the ion data is not taken into account. This option makes it possible to correlate ion distributions with given helical or backbone parameter states.

7) Root mean square fluctuations:

Setting **rmsf**=t. uses the combination of the helical ion parameters and an average helical axis to map the ions into Cartesian space and then calculates their average position (*lis file_cen.pdb*) and their root mean square fluctuation values (*lis file.rmsf*). A single pass rmsf algorithm to make this calculation possible with a single read of the trajectory file. This option is generally used for solute atoms and not for solvent molecules or ions.

8) **circ**:

Should be set to .t. when minicircles are analyzed. In this case, seq can be used to search sequences that overlap the minicircle junction (base pairs N and 1).

9) **istep**:

Outputs the average ion population (in a **lis.stp** file) every **istep** snapshots. This is useful for calculating the convergence of an ion population in a volume specified by a combination of dlow/dhig/rlow/rhig/allow/ahig values. **istep** is not compatible with .cub input data.

10) **seqin**:

Inputting the first strand base sequence of the oligomer analyzed via **seqin** is necessary when reading .cub or .pts data.

11) **Series output**:

If **series** is true, the snapshot number and the number of ions counted for the snapshot is output in a file (**lis.cser**). This option is generally useful when a limited zone of space is analyzed.

13) Taking **solute** volume into account:

By reading a pdb file containing the nucleic acid fragment (and nothing else), it is possible to take into account the space occupied by the van der Waals volume of the nucleic acid and to eliminate this volume when calculating molarities. This gives a better view of the balance of ion densities between the major and minor grooves, where the space occupied by the nucleic acid is different. Note that the solute conformation should be consistent with the axfrm data describing its helical axis. We find that for a trajectory analysis it is generally best to use the average structure obtained after superposing the snapshots along the trajectory (rather than using the snapshot with the smallest rmsd to this average structure).

Other input data options

There are two alternatives to reading .cdi data from Curves+

1) .pts input

This is a simplified input where a single file contains an unlimited number of snapshots. Each snapshot starts with the number of ions to be read (format i10). this is followed by the Cartesian coordinates of the ions as a list of x,y,z values (format 3f8.3). Note that this input does not allow more the one type of ion to be read and so **type** should be left at its default value "*".

2) .cub input

This allows a cubic density matrix to be read in Gaussian cube format. The format begins with two title lines (any text), then a line specifying the x, y, z coordinates of the lower corner of the cubic space (format 4x, 3f10.3) and three lines specifying the number of points along each axis (respectively, idx, iyd, izd) and a unit vector defining the axis (format i4, 3f10.3). The density matrix *his* is then read in free format:

```
do i=1,idx
do j=1,iyd
read(2,*) (his(I,j,k), k=1,izd)
enddo
enddo
```

Output files:

CANION analyzes ion/atom distributions in 1D, 2D and 3D

1D) histograms of distributions with respect to **d**, **r**, and **a** (in files .d, .r and .a). In the **r** histogram, each bar is divided by **r** to correct for the increasing radial volume being analyzed. Expressed in molarity.

2D) histograms of the distributions with respect to **dr**, **da** and **ra** (in files .dr, .da and .ra). The **dr** distribution is corrected for increasing radial volume as described above. The **ra** distribution is plotted in polar coordinates. Expressed in molarity.

3D) an ion/atom density matrix is generated with respect to the average (or standard) helical axis using **grid** to set the spacing and **shell** to determine the size of the matrix. Note that you can use any axis for this analysis (as long as it has the correct length). For example, you can analyze a minicircle using a straight axis. The output format is a Gaussian density matrix, expressed in units of molarity.

.ria) Accumulated ion population as a function of radius **r**.

.stp) Average ion population at chosen time steps (see namelist option **istep**).

.cser) number of ions in chosen volume at each snapshot (see namelist option **series**).

Sample I/P file for CANION:

```
rm Ant_KG.*
/Users/rl/Code/util/canion <<!
  &inp dat=Antion4.cdi,axfrm=caver,solute=aver,
  lis=Ant_KG,seq=G,type=K+, &end
!
```

- this run analyzes potassium ions around GC base pairs (using the ion position data in Antion4.cdi and the average DNA conformation in aver.pdb, with corresponding helical axis frames in caver.afr)

Appendix: MatLab scripts for analysis

MatLab offers a convenient way to display output from CANION. Here are some functions for plotting 1D and 2D data. For 3D data, use molecular graphic programs such as CHIMERA (<http://www.cgl.ucsf.edu/chimera/>) to read the .cub files output by Canion.

1D plot: using Matlab command `csing('file.r',[slen])`

- where the file extension can be .r, .d or .a. The length of the nucleic acid fragment is assumed to be 18, but it can be reset with the optional variable `slen`

```
function x=csing(fir,slen)
a1=33;
a2=147;
pr=10.25;
siz=18;
    if nargin==2
        siz=slen;
    end
lw=2.0;
head=['CANION 1D: ',fir];
scr=get(0,'ScreenSize');
figure('Position',[scr(3)/2 scr(4)/2 scr(3)/2 scr(4)/2],'Name',head)
fr=load(fir);
flim=size(fr,1);
plot(fr(1:flim,1),fr(1:flim,2),'k','LineWidth',2)
set(gca,'FontSize',18);
ylabel('Molarity','FontSize',18);
limy=ylim;
    hold on;
    if size(strfind(fir,'.r'),1)==1
        ylim(limy);
        plot([pr;pr],limy,'k','LineWidth',lw);
        xlabel('R (\AA)','Interpreter','Latex','FontSize',18);
    elseif size(strfind(fir,'.a'),1)==1
        ylim(limy);
        xlim([1,360]);
        plot([a1;a1],limy,'-k','LineWidth',lw);
        plot([a2;a2],limy,'-k','LineWidth',lw);
        xlabel('A (degrees)','FontSize',18);
    elseif size(strfind(fir,'.d'),1)==1
        ylim(limy);
        xlim([1,siz]);
        xlabel('D (bp)','FontSize',18);
    end
end
```

Comparing two 1D plots: using Matlab command `csing('file1.r','file2.r', [slen])`. Note that the y-axis limit is set by the maximum y-value in file1.r.

```
function x=csing2(fir,sec,slen)
a1=33;
a2=147;
pr=10.25;
siz=18;
    if nargin==3
        siz=slen;
    end
lw=2.0;
head=['CANION 1D: ',fir];
scr=get(0,'ScreenSize');
figure('Position',[scr(3)/2 scr(4)/2 scr(3)/2 scr(4)/2],'Name',head)
fr=load(fir);
sc=load(sec);
plot(fr(:,1),fr(:,2),'k','LineWidth',2)
limy=ylim;
ylabel('Molarity','FontSize',18);
hold on;
plot(sc(:,1),sc(:,2)-min(sc(:,2)),'--k','LineWidth',2)
    if size(strfind(fir,'.r'),1)==1
        ylim(limy);
        plot([pr;pr],limy,'--k','LineWidth',lw);
        xlabel('R (\AA)','Interpreter','Latex','FontSize',18);
    elseif size(strfind(fir,'.a'),1)==1
        ylim(limy);
        xlim([1,360]);
        plot([a1;a1],limy,'k','LineWidth',lw);
        plot([a2;a2],limy,'k','LineWidth',lw);
        xlabel('A (degrees)','FontSize',18);
    elseif size(strfind(fir,'.d'),1)==1
        ylim(limy);
        xlim([1,siz])
        xlabel('D (bp)','FontSize',18);
    end
end
```


2D plots: using Matlab command `csurf('file',{rlim})` - produces three plots for the planes *DA*, *DR* and *RA* from the corresponding data files. The optional variable *rlim* can be used to zoom in on a smaller range of *R* in the *DR* and *DA* plots. *rlim* must be \leq *rhig*, where *rhig* was the upper limit of the radius used in the corresponding Canion calculations.

```
function x=csurf(fir,rlim)
head=['CANION 2D: ',fir];
scr=get(0,'ScreenSize');
fda=load([fir,'.da']);
fdr=load([fir,'.dr']);
fra=load([fir,'.ra']);
kpd=6;
kpr=2;
kpa=0.2;
lw=2.0;
pr=10.25;
a1=33; ar1=a1*pi/180;
a2=147; ar2=a2*pi/180;
dx=size(fda,1);
ax=size(fda,2);
rx=size(fra,1);
dlow=1+0.5/kpd;
dhig=dlow+(dx-1)/kpd;
rlow=0.5/kpr;
rhig=rlow+(rx-1)/kpr;
    if nargin==2
        rhig=rlim;
        rx=round(rlow+kpr*rhig+1);
        fdr=fdr(:,1:rx);
        fra=fra(1:rx,:);
    end
alow=0.5/kpa;
ahig=alow+(ax-1)/kpa;
d=linspace(dlow,dhig,dx);
r=linspace(rlow,rhig,rx);
a=linspace(alow,ahig,ax);

figure('Position',[scr(3)/9, scr(4)/2 scr(3)/2.5
scr(4)/2],'Name',head);
contourf(a,d,fda);
set(gca,'FontSize',18);
colorbar('FontSize',18);
axis square;
title('DA','FontSize',18);
ylabel('D (base pair steps)','FontSize',18);
xlabel('A (degrees)','Interpreter','Latex','FontSize',18);
limy=ylim;
    hold on;
    ylim(limy);
    xlim([1,360]);
    plot([a1;a1],limy,'-w','LineWidth',lw);
    plot([a2;a2],limy,'-w','LineWidth',lw);
figure('Position',[2*scr(3)/9, scr(4)/2 scr(3)/2.5
scr(4)/2],'Name',head);
contourf(r,d,fdr);
set(gca,'FontSize',18);
colorbar('FontSize',18);
axis square;
title('DR','FontSize',18);
ylabel('D (bp)','FontSize',18);
```

```

xlabel('R (\AA)', 'Interpreter', 'Latex', 'FontSize', 18);
limy=ylim;
    hold on
    ylim(limy);
    plot([pr;pr], limy, '-w', 'LineWidth', lw);
figure('Position', [3*scr(3)/9, scr(4)/2 scr(3)/2.5
scr(4)/2], 'Name', head);
[theta,r] = meshgrid(linspace(-pi,pi,ax),linspace(rlow,rhig,rx));
contourf(r.*cos(theta),r.*sin(theta),fra);
set(gca, 'FontSize', 18);
ylabel('Y (\AA)', 'Interpreter', 'Latex', 'FontSize', 18);
xlabel('X (\AA)', 'Interpreter', 'Latex', 'FontSize', 18);
    hold on
    plot([-cos(ar1)*pr;0;-cos(ar2)*pr],[-sin(ar1)*pr;0;-sin(ar2)*pr],
        '-w', 'LineWidth', lw);
    plot([0;0],[0;pr], '-w', 'LineWidth', lw);
    ang=0:0.01:2*pi;
    xp=pr*cos(ang);
    yp=pr*sin(ang);
    plot(xp,yp, '-w', 'LineWidth', lw);
    axis square;
    axis tight;
    colorbar('FontSize', 18);
    title('RA', 'FontSize', 18);

```